


Ceramic Evaporation Boat Set Alumina Crucible For Laboratory Use

Item Number: KME09



Introduction

It can be used for vapor deposition of various metals and alloys. Most metals can be evaporated completely without loss. Evaporation baskets are reusable.¹

[Learn More](#)

Material	ceramic crucible / tungsten wire
Ceramic Crucible Size	0.25ml / 0.5ml / 1.0ml / 2.0ml / 3.0ml
Tungsten Wire Size	3*0.8 wire / 3*1.0 wire / 3*0.8 wire