

## Glass Tissue Grinder Homogenizer For Lab Use

Item Number: KTG-1



### Introduction

The tiny gaps and slight bumps between the glass tube and the plunger are engaged, and a grinding effect is produced when the plunger is rotated.

[Learn More](#)

Model	Total length (mm)	Frosted length (mm)	Outer diameter of tube (mm)	Inner diameter of mouth (mm)
0.5ml	133	34	6	16
1ml	147	47	8	17
2ml	155	53.5	10	20
3ml	160	61.5	10	20
5ml	195	72	11	24
10ml	225	74	13	26
15ml	245	83	16	28
20ml	255	92	17	31.5
25ml	270	110	18	37
30ml	275	112	18	34
50ml	295	112	21	38