

Laboratory Formalin Specimen Jar Sealed Sample Bottle

Item Number: KTG-3



Introduction

Specimen bottles are used to collect and preserve animal and plant specimens.

[Learn More](#)

Specifications	Overall height	Bottle height	Bottle diameter	Bottom diameter	Bottle mouth outer diameter
45*65mm	105mm	78mm	45	59mm	44
45*90mm	132mm	106mm	45	59mm	54
45*120mm	168mm	136mm	45	58mm	54
45*150mm	189mm	165mm	45	60mm	56
60*90mm	135mm	105mm	60	76mm	68
60*120mm	170mm	135mm	60	75mm	68
60*150mm	200mm	166mm	60	75mm	69
60*180mm	228mm	196mm	60	74mm	69
75*90mm	140mm	104mm	75	95mm	85
75*120mm	172mm	135mm	75	92mm	87
75*150mm	200mm	166mm	75	93mm	87
75*180mm	230mm	198mm	75	94mm	84
90*90mm	146mm	105mm	90	109mm	98
90*120mm	183mm	136mm	90	109mm	98
90*150mm	205mm	165mm	90	111mm	101
90*180mm	233mm	198mm	90	109mm	98