

Vibrating Disc Mill Small Laboratory Grinding Machine

Item Number: KT-JZ



Introduction

Discover the versatile Vibrating Disc Mill for efficient laboratory grinding. Ideal for geology, metallurgy, biology, and more. Explore now!

[Learn More](#)

| Serial Number | Model | Feed (mm) | Discharge (mesh) | Loading Amount (g) | Quantity | Voltage (V) | Pot Voltage (kW) | Dimensions (LxWxH, mm) | Weight (kg) |
|---------------|-------------|-----------|------------------|--------------------|----------|-------------|------------------|------------------------|-------------|
| 1 | JZGJ-50-1 | | 80-300 | 50 | 1 | 220 | 0.18 | 300x300x400 | 45 |
| 2 | JZGJ-50-1D | | 80-300 | 80-300 | 1 | 220 | 0.18 | 300x300x400 | 45 |
| 3 | JZGJ-80-1 | | 80-300 | 80 | 1 | 220 | 0.25 | 380x380x600 | 65 |
| 4 | JZGJ-80-1D | | 80-300 | 80 | 1 | 220 | 0.25 | 380x380x600 | 65 |
| 5 | JZGJ-80-2 | | 80-300 | 80x2 | 2 | 220 | 0.25 | 380x380x600 | 80 |
| 6 | JZGJ-80-2D | | 80-300 | 80x2 | 2 | 220 | 0.25 | 380x380x600 | 80 |
| 7 | JZGJ-100-1 | | 80-300 | 100 | 1 | 380 | 1.1 | 540x520x780 | 130 |
| 8 | JZGJ-100-1D | | 80-300 | 100 | 1 | 380 | 1.1 | 540x520x780 | 130 |
| 9 | JZGJ-100-2 | | 80-300 | 100x2 | 2 | 380 | 1.1 | 540x520x780 | 160 |
| 10 | JZGJ-100-2D | | 80-300 | 100x2 | 2 | 380 | 1.1 | 540x520x780 | 160 |
| 11 | JZGJ-100-3 | | 80-300 | 100x3 | 3 | 380 | 1.5 | 550x600x850 | 200 |
| 12 | JZGJ-100-3D | | 80-300 | 100x3 | 3 | 380 | 1.5 | 550x600x850 | 200 |
| 13 | JZGJ-100-4 | | 80-300 | 100x4 | 4 | 380 | 2.2 | 600x550x850 | 210 |
| 14 | JZGJ-100-4D | | 80-300 | 100x4 | 4 | 380 | 2.2 | 600x550x850 | 210 |
| 15 | JZGJ-100-6 | | 80-300 | 100x6 | 6 | 380 | 2.2 | 700x650x900 | 250 |
| 16 | JZGJ-100-6D | | 80-300 | 100x6 | 6 | 380 | 2.2 | 700x650x900 | 250 |
| 17 | JZGJ-100-7 | | 80-300 | 100x7 | 7 | 380 | 2.2 | 700x650x900 | 260 |
| 18 | JZGJ-100-7D | | 80-300 | 100x7 | 7 | 380 | 2.2 | 700x650x900 | 260 |
| 19 | JZGJ-200-1 | | 80-300 | 200 | 1 | 380 | 1.5 | 540x520x780 | 180 |
| 20 | JZGJ-200-1D | | 80-300 | 200 | 1 | 380 | 1.5 | 540x520x780 | 180 |
| 21 | JZGJ-200-2 | | 80-300 | 200x2 | 2 | 380 | 1.5 | 550x600x850 | 210 |
| 22 | JZGJ-200-2D | | 80-300 | 200x2 | 2 | 380 | 1.5 | 550x600x850 | 210 |
| 23 | JZGJ-200-3 | | 80-300 | 200x3 | 3 | 380 | 2.2 | 700x650x900 | 220 |
| 24 | JZGJ-200-3D | | 80-300 | 200x3 | 3 | 380 | 2.2 | 700x650x900 | 220 |
| 25 | JZGJ-400-1 | | 80-300 | 400 | 1 | 380 | 1.5 | 600x550x850 | 190 |

| Serial Number | Model | Feed (mm) | Discharge (mesh) | Loading Amount (g) | Quantity | Voltage (V) | Pot Voltage (kW) | Dimensions (LxWxH, mm) | Weight (kg) |
|---------------|--------------|-----------|------------------|--------------------|----------|-------------|------------------|------------------------|-------------|
| 26 | JZGJ-400-1D | | 80-300 | 400 | 1 | 380 | 1.5 | 600x550x850 | 190 |
| 27 | JZGJ-1000-1D | | 80-300 | 1000 | 1 | 380 | 2.2 | 700x650x900 | 260 |
| 28 | JZGJ-100-1DB | | 80-300 | 100 | 1 | 380 | 1.1 | 700x520x800 | 170 |
| 29 | JZGJ-100-2DB | | 80-300 | 100x2 | 2 | 380 | 1.1 | 700x520x800 | 180 |
| 30 | JZGJ-100-3DB | | 80-300 | 100x3 | 3 | 380 | 2.2 | 1000x580x900 | 210 |
| 31 | JZGJ-100-4DB | | 80-300 | 100x4 | 4 | 380 | 2.2 | 1000x580x900 | 240 |
| 32 | JZGJ-200-1DB | | 80-300 | 200 | 1 | 380 | 1.5 | 700x520x800 | 190 |