

1-5L Jacketed Glass Reactor Vessel Lab Glass Reactor

Item Number: KR-1



Introduction

Discover the perfect solution for your pharmaceutical, chemical, or biological products with our 1-5L jacket glass reactor system. Custom options available.

[Learn More](#)

Model	KGR-1	KGR-2	KGR-3	KGR-5
Reactor Volume	1L	2L	3L	5L
Jacket Volume	0.3L	0.6L	0.9L	1.5L
Funnel size	0.2L	0.2L	0.2L	0.5L
Vacuum Sealing	Teflon+Viton Sealing			
Speed Control	Knop Stepless Speed Regulation			
Stirring Speed	0-1300rpm			
Working temperature	-80~250°C			
Vacuum	0.098Mpa			
Safety Protection	Over-current/ ground fault/over-temperature protection			
Working Environment	5-35°C			
Power Supply	110-440V, 50/60Hz			
Overall Dimensions	480*420*1110 mm	480*420*1110 mm	480*420*1110 mm	591*520*1594mm
Net weight	19kg	20kg	21kg	34kg