

# 100L Chilling Circulator Cooling Water Circulator For Low Temperature Constant Temperature Reaction Bath Water Bath Cooling

Item Number: KCP-100



## Introduction

Get reliable and efficient chilling power for your lab or industrial needs with KinTek KCP chilling circulator. With max. -120°C temperature and built-in circulating pump.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP100-20	100L	-23-5°C	9738-2626W	70L/min 6-8m	+/-2°C	1/2"	1100*900*1100
KCP100-30		-32-5°C	9820-2820W				
KCP100-40		-40-5°C	10860-2830W				
KCP100-60		-63-5°C	11380-2100W				1200*1100*1350
KCP100-80		-83-5°C	11670-1800W				
KCP100-120		-120-5°C	11880-1200W				