

10L Chilling Circulator Cooling Water Bath Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-10



Introduction

Get the KinTek KCP 10L Chilling Circulator for your lab needs. With a stable and quiet chilling power of up to -120°C, it also works as a one chilling bath for versatile applications.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP10-10	10L	-12-5°C	1825-820W	35L/min 4-6m	+/-2°C	1/2"	580*450*730
KCP10-20		-23-5°C	2100-780W				
KCP10-25		-26-5°C	2170-780W				
KCP10-30		-32-5°C	2436-820W				720*650*900
KCP10-40		-40-5°C	2548-980W				
KCP10-60		-63-5°C	2732-310W				
KCP10-80		-83-5°C	2900-230W				
KCP10-120		-120-5°C	2982-200W				940*710*900