

10L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBH-10



Introduction

Experience efficient lab performance with KinTek KCBH 10L Heating Chilling Circulator. Its all-in-one design offers reliable heating, chilling, and circulating functions for industrial and lab use.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)
KCBH10-10	10L	-10-200°C	1865-860W	7Kw	20L/min 20m	+/-2°C	540*440*800
KCBH10-20		-20-200°C	2140-359W				
KCBH10-30		-30-200°C	1415-490W				
KCBH10-40		-40-200°C	1460-501W	7.5Kw			630*550*1400
KCBH10-60		-60-200°C	1950-270W				
KCBH10-80		-80-200°C	1978-240W				