

2-5L Rotary Evaporator

Item Number: KE-5



Introduction

Efficiently remove low boiling solvents with the KT 2-5L Rotary Evaporator. Perfect for chemical labs in the pharmaceutical, chemical, and biological industries.

[Learn More](#)

| Model | KRE-2000 | KRE3000 | KRE5000 |
|----------------------|--|--------------|--------------|
| Evaporating Flask | 2L | 3L | 5L |
| Recycling Flask | 1L | 2L | 3L |
| Condenser | Vertical/ Diagonal | | |
| Vacuum Degree | ≤2kPa/h | | |
| Vacuum Sealing | Teflon+Viton Sealing | | |
| Speed Control | Knop Stepless Speed Regulation | | |
| Rotating Speed | 5-200rpm | | |
| Rotating Power | 40W | | |
| Heating Power | 1.4KW | 1.5KW | 1.8Kw |
| Heating Range | 0-99°C | | |
| Temperature Accuracy | ±1°C | | |
| Lifting Method | Electric | | |
| Lifting Distance | 150mm | 180mm | 180mm |
| Anti-splash Device | PP Transparent Cover (Optional) | | |
| Safety Protection | Over-current/ ground fault/over-temperature protection | | |
| Power supply | 110-240V, 50/60Hz | | |
| Working Environment | 5-35°C | | |
| Dimensions | 680*320*760mm | 550*420*78mm | 550*420*78mm |
| Weight | 25Kg | 35kg | 36kg |