

## 20L Chiller Water Bath Cooling Circulator Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-20



### Introduction

KinTek KCP chilling circulator is a versatile and reliable equipment that supplies constant chilling power with circulating fluids. It can work as a one chilling bath and reach a max. Chilling temperature of -120°C.

[Learn More](#)

| Model     | Fluid Volume | Working Temp. | Cooling Capacity | Pump Flowrate | Temp. Stability | Hose Port | Dimensions (mm) |
|-----------|--------------|---------------|------------------|---------------|-----------------|-----------|-----------------|
| KCP20-20  | 20L          | -23-5°C       | 3126-950W        | 35L/min 4-6m  | +/-2°C          | 1/2"      | 630*52*1000     |
| KCP20-30  |              | -32-5°C       | 3237-1023W       |               |                 |           |                 |
| KCP20-40  |              | -40-5°C       | 3635-980W        |               |                 |           |                 |
| KCP20-60  |              | -63-5°C       | 4210-510W        |               |                 |           | 940*710*1010    |
| KCP20-80  |              | -83-5°C       | 4530-360W        |               |                 |           |                 |
| KCP20-120 |              | -120-5°C      | 4889-270W        |               |                 |           |                 |