

20L Chilling Circulator

Item Number: KCP-20



Introduction

KinTek KCP chilling circulator is a versatile and reliable equipment that supplies constant chilling power with circulating fluids. It can work as a one chilling bath and reach a max. Chilling temperature of -120°C.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP20-20	20L	-23-5°C	3126-950W	35L/min 4-6m	+/-2°C	1/2"	630*52*1000
KCP20-30		-32-5°C	3237-1023W				
KCP20-40		-40-5°C	3635-980W				
KCP20-60		-63-5°C	4210-510W				940*710*1010
KCP20-80		-83-5°C	4530-360W				
KCP20-120		-120-5°C	4889-270W				