

20L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBH-20



Introduction

Maximize lab productivity with KinTek KCBH 20L Heating Chilling Circulator. Its all-in-one design offers reliable heating, chilling, and circulating functions for industrial and lab use.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)
KCBH20-10	20L	-10-200°C	3126-950W	7.2Kw	42L/min 28m	+/-2°C	630*550*1100
KCBH20-20		-20-200°C	3163-980W				
KCBH20-30		-30-200°C	3277-1061W				
KCBH20-40		-40-200°C	3675-998W	7.8Kw			
KCBH20-60		-60-200°C	4160-550W				
KCBH20-80		-80-200°C	4576-420W				