

20L Rotary Evaporator

Item Number: KE-20



Introduction

Efficiently separate "low boiling" solvents with the 20L Rotary Evaporator, ideal for chemical labs in pharmaceutical and other industries. Guarantees working performance with selected materials and advanced safety features.

[Learn More](#)

Model	KRE-2020Z
Evaporating Flask	20L
Flask Size	ø125mm
Receiving Flask	10L
Condenser	Vertical/Three Flow
Vacuum Degree	≤2kPa/h
Vacuum Sealing	Teflon+Viton Sealing
Speed Control	Knop Stepless Speed Regulation
Rotating Speed	10-120rpm
Rotating Power	125W
Heating Power	4KW
Heating Range	0-99°C
Temperature Accuracy	±2°C
Lifting Method	Electric
Lifting Distance	200mm
Anti-splash Device	PP Transparent Cover (Optional)
Safety Protection	Over-current/ ground fault/over-temperature protection
Power supply	110-240V, 50/60Hz
Working Environment	5-35°C
Dimensions	850*520*160mm
Weight	90Kg