

23L Laboratory Vacuum Drying Oven

Item Number: KT-VDO



Introduction

Kintek Intelligent Vacuum Drying Oven for labs: precise, stable, low-temp drying. Ideal for heat-sensitive materials. Get a quote now!

[Learn More](#)

Parameter	Specification
Model	KT-DZK-6025B
Display Mode	LCD Screen
Temp. Control Range	RT (Room Temp) +10°C ~ 250°C
Temp. Resolution	0.1°C
Vacuum Level	-0.1 MPa
Timer Range	1 min ~ 999 hr 59 min
Heating Power	400W (Total Power: 800W)
Max Continuous Run	7 Days (at 230°C)
Chamber Material	201 Stainless Steel (304 SS Optional)
Chamber Dimensions	280 (D) × 280 (W) × 300 (H) mm
Capacity	23 L
Shelves	2 Standard Shelves (Load: 15kg each)
Power Supply	AC 220V / 50Hz
Net Weight	46 Kg