

# 30L Chiller Water Bath Cooling Circulator Low Temperature Constant Temperature Reaction Bath


Item Number: KCP-30



## Introduction

Keep your lab cool with the KinTek KCP chilling circulator - ideal for constant chilling power and adaptable to meet all your working needs.

[Learn More](#)



Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP30-20	30L	-23-5°C	3780-1150W	35L/min 4-6m	+/-2°C	1/2"	940*710*1010
KCP30-30		-32-5°C	3960-1200W				
KCP30-40		-40-5°C	4080-1270W				
KCP30-60		-63-5°C	4330-680W				1030*810*1020
KCP30-80		-83-5°C	4679-422W				
KCP30-120		-120-5°C	4968-330W				