

# 30L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBH-30



## Introduction

Get versatile lab performance with KinTek KCBH 30L Heating Chilling Circulator. With max. heating temp of 200°C and max. chilling temp of -80°C, it's perfect for industrial needs.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)
KCBH30-10	30L	-10-200°C	3780-1150W	10.2Kw	42L/min 28m	+/-2°C	630*550*1100
KCBH30-20		-20-200°C	3960-1200W				
KCBH30-30		-30-200°C	7803-3746W				
KCBH30-40		-40-200°C	9931-4190W	12.5Kw			
KCBH30-60		-60-200°C	10438-3219W				
KCBH30-80		-80-200°C	10578-3087W				