

40L Chiller Water Bath Cooling Circulator Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-40



Introduction

Get efficient and reliable chilling power with KinTek KCP circulating chiller. With a max. temp of -120°C, it's an ideal equipment for different working circumstances.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP40-20	40L	-23-5°C	3780-1150W	70L/min 6-8m	+/-2°C	1/2"	1030*860*1030
KCP40-30		-32-5°C	3960-1200W				
KCP40-40		-40-5°C	4080-1270W				
KCP40-60		-63-5°C	4330-680W				
KCP40-80		-83-5°C	4679-422W				
KCP40-120		-120-5°C	4968-330W				