

50L Chiller Water Bath Cooling Circulator Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-50



Introduction

KinTek KCP 50L chilling circulator is a reliable and efficient equipment for supplying constant chilling power with circulating fluids in various working circumstances.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP50-20	50L	-23-5°C	8520-2170W	70L/min 6-8m	+/-2°C	1/2"	760*610*1030
KCP50-30		-32-5°C	8930-2280W				
KCP50-40		-40-5°C	10256-23300W				
KCP50-60		-63-5°C	10570-1600W				1030*860*1030
KCP50-80		-83-5°C	11200-1480W				
KCP50-120		-120-5°C	11290-1260W				