

50L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBH-50



Introduction

Experience versatile heating, chilling, and circulating capabilities with our KinTek KCBH 50L Heating Chilling Circulator. Ideal for labs and industrial settings, with efficient and reliable performance.

[Learn More](#)

| Model | Fluid Volume | Working Temp. | Cooling Capacity | Overall Power | Flowrate | Temp. Stability | Dimensions (mm) |
|-----------|--------------|---------------|------------------|---------------|-------------|-----------------|-----------------|
| KCBH50-20 | 50L | -20-200°C | 11156-5341W | 12 | 42L/min 28m | +/-2°C | 720*650*1150 |
| KCBH50-30 | | -30-200°C | 13695-9575W | | | | |
| KCBH50-40 | | -40-200°C | 11743-7351W | | | | |
| KCBH50-60 | | -60-200°C | 18437-7601W | 15.5Kw | | | 880*880*1250 |
| KCBH50-80 | | -80-200°C | 18657-7760W | | | | |