

# 50L Heating Chilling Circulator

Item Number: KCBH-50



## Introduction

Experience versatile heating, chilling, and circulating capabilities with our KinTek KCBH 50L Heating Chilling Circulator. Ideal for labs and industrial settings, with efficient and reliable performance.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)
KCBH50-20	50L	-20-200°C	11156-5341W	12	42L/min 28m	+/-2°C	720*650*1150
KCBH50-30		-30-200°C	13695-9575W				
KCBH50-40		-40-200°C	11743-7351W				
KCBH50-60		-60-200°C	18437-7601W	15.5Kw			880*880*1250
KCBH50-80		-80-200°C	18657-7760W				