

5L Chilling Circulator Cooling Water Bath Circulator For Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-5



Introduction

Maximize lab efficiency with the KinTek KCP 5L Chilling Circulator. Versatile and reliable, it provides constant chilling power up to -120°C.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP5-10	5L	-12-5°C	1213-290W	35L/min 4-6m	+/-2°C	1/2"	520*350*720
KCP5-20		-23-5°C	1248-319W				
KCP5-25		-26-5°C	1314-447W				
KCP5-30		-32-5°C	1375-450W				
KCP5-40		-40-5°C	1420-462W				430*350*960
KCP5-60		-63-5°C	1910-230W				
KCP5-80		-83-5°C	1938-200W				
KCP5-120		-120-5°C	2635-180W				