

5L Chilling Circulator

Item Number: KCP-5



Introduction

Maximize lab efficiency with the KinTek KCP 5L Chilling Circulator. Versatile and reliable, it provides constant chilling power up to -120°C.

[Learn More](#)

| Model | Fluid Volume | Working Temp. | Cooling Capacity | Flowrate | Temp. Stability | Hose Port | Dimensions (mm) |
|----------|--------------|---------------|------------------|--------------|-----------------|-----------|-----------------|
| KCP5-10 | 5L | -12-5°C | 1213-290W | 35L/min 4-6m | +/-2°C | 1/2" | 520*350*720 |
| KCP5-20 | | -23-5°C | 1248-319W | | | | |
| KCP5-25 | | -26-5°C | 1314-447W | | | | |
| KCP5-30 | | -32-5°C | 1375-450W | | | | |
| KCP5-40 | | -40-5°C | 1420-462W | | | | 430*350*960 |
| KCP5-60 | | -63-5°C | 1910-230W | | | | |
| KCP5-80 | | -83-5°C | 1938-200W | | | | |
| KCP5-120 | | -120-5°C | 2635-180W | | | | |