

5L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBH-5



Introduction

KinTek KCBH 5L Heating Chilling Circulator - Ideal for labs and industrial conditions with multi-functional design and reliable performance.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)	
KCBH5-10	5L	-10-200°C	1243-330W	3.5Kw	20L/min 20m	+/-2°C	540*440*800	
KCBH5-20		-20-200°C	1288-359W					
KCBH5-30		-30-200°C	1415-490W	4Kw				
KCBH5-40		-40-200°C	1460-501W					
KCBH5-60		-60-200°C	1950-270W	7Kw				630*550*1400
KCBH5-80		-80-200°C	1978-240W					