

80L Chilling Circulator Cooling Water Circulator For Water Bath Cooling And Low Temperature Constant Temperature Reaction Bath

Item Number: KCP-80



Introduction

Efficient and Reliable 80L Chilling Circulator with a max temp of -120°C. Ideal for labs and industrial use, also works as a one chilling bath.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Pump Flowrate	Temp. Stability	Hose Port	Dimensions (mm)
KCP80-20	80L	-23-5°C	9738-2626W	70L/min 6-8m	+/-2°C	1/2"	1000*800*1100
KCP80-30		-32-5°C	9820-2820W				
KCP80-40		-40-5°C	10860-2830W				
KCP80-60		-63-5°C	11380-2100W				1100*830*1250
KCP80-80		-83-5°C	11670-1800W				
KCP80-120		-120-5°C	11880-1200W				