

# 80L Heating Chilling Circulator Cooling Water Bath Circulator For High And Low Temperature Constant Temperature Reaction

Item Number: KCBP-80



## Introduction

Get all-in-one heating, chilling, and circulating capabilities with our KinTek KCBH 80L Heating Chilling Circulator. High efficiency, reliable performance for labs and industrial applications.

[Learn More](#)

Model	Fluid Volume	Working Temp.	Cooling Capacity	Overall Power	Flowrate	Temp. Stability	Dimensions (mm)
KCBH100-20		-20-200°C	18638-6966W				
KCBH100-30	100L	-30-200°C	19904-9557W	13.5	42L/min 28m	+/-2°C	880*880*1250
KCBH100-40		-40-200°C	25334-10698W				