

Laboratory Manual Slicer

Item Number: PT10



Introduction

Manual microtome is a high-precision cutting device designed for laboratories, industry and medical fields. It is suitable for the preparation of thin slices of various materials such as paraffin samples, biological tissues, battery materials, food, etc.

[Learn More](#)

Instrument model	PT10
Mold material	3Cr13 (quenched stainless steel)
Slice size	Φ10□Φ12□Φ14□Φ15□Φ16□Φ18□Φ19(d)
Overall dimensions	120*200*280mm(L*W*H)
Equipment weight	10kg
Machine size diagram	