

Anti-Cracking Press Mold For Lab Use

Item Number: PML



Introduction

The anti-cracking press mold is a specialized equipment designed for molding various shapes and sizes of film using high pressure and electric heating.

[Learn More](#)

Instrument model	PMT
Sample shape	
Die heating	Alloy tool steel :Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ6□ Φ8□ Φ10□ Φ15□ Φ20mm (M)
Cavity depth	40mm (N)
Dimensions	Φ98*120mm(L*H)
Weight	5Kg
Diagram of hydraulic powder press size	