

Assemble Lab Cylindrical Press Mold

Item Number: PMAC



Introduction

Get reliable and precise molding with Assemble Lab Cylindrical Press Mold. Perfect for ultra-fine powder or delicate samples, widely used in material research and development.

[Learn More](#)

Model	PMAS-A	PMAS-B	PMAS-C	PMAS-D	PMAS-E
Material	Cr12MoV				
Indenter hardness	HRC60-HRC62				
Sample size	Φ3[]Φ4[]Φ5[]Φ6[]Φ8[]Φ10mm (M)	Φ12[]Φ13[]Φ15[]Φ18[]Φ20mm (M)	Φ30[]Φ40mm (M)	Φ50[]Φ60mm (M)	Φ70[]Φ80 mm (M)
Cavity depth	30mm (N)	40 mm (N)	50 mm (N)	55 mm (N)	60 mm (N)
External dimensions	Φ43*93mm(L*H)	Φ53*120mm(L*H)	Φ73*133mm[]Φ95*133mm(L*H)	Φ115*150mm[]Φ127*150mm (L*H)	Φ153*180mm[]Φ180*180mm(L*H)
Weight(Kg)	0.75Kg	1.2Kg	3.8Kg[]6.3Kg	14Kg[]20Kg	30Kg[]40Kg

Other sizes can be customized