

Ball Press Mold For Lab

Item Number: PMQ



Introduction

Explore versatile Hydraulic Hot Press molds for precise compression molding. Ideal for creating various shapes and sizes with uniform stability.

[Learn More](#)

Instrument model	PMQ	
Sample shape		
Die heating	Alloy tool steel :Cr12MoV	
Indenter hardness	HRC60-HRC62	
Sample size	Φ6□Φ8□Φ10□Φ15□Φ20mm (M)	Φ30□Φ40□Φ50 mm (M)
Cavity depth	40mm (N)	60 mm (N)
Dimensions	Φ53*120mm (L*H)	Φ88*150 mm (L*H)
Weight	1.4kg	5.8kg
Diagram of hydraulic powder press size		