

Benchtop High Speed Small Lab Centrifuge For Serum Separation And Low Temperature Dehydration

Item Number: KT-SCS



Introduction

Kintek Benchtop High-Speed Centrifuge: Efficient, stable, precise. Ideal for clinical, research labs. Explore features & specs.

[Learn More](#)

Parameter	Specification
Display	LCD Display
Max. Speed	18,500 rpm
Speed Control Accuracy	±10 rpm
Max. RCF	22,575 ×g
Max. Capacity	6 × 100 mL
Chamber Diameter	Φ280 mm
Timer Range	1s to 99min 59s
Noise Level	≤55 dB
Power Supply	AC 220V±22V 50Hz
Power Consumption	400W
Net Weight	28 kg
Dimensions (L×W×H)	360 × 470 × 320 mm