

# Benchtop Water Circulating Vacuum Pump

Item Number: KV-2



## Introduction

Need a water circulating vacuum pump for your lab or small-scale industry? Our Benchtop Water Circulating Vacuum Pump is perfect for evaporation, distillation, crystallization, and more.

[Learn More](#)

### Detailed material

Ejector	check valve	Connecting pipe	three links	pump body	impeller	Air extraction nozzle
PP	Copper+PP	NR	PP	PP	PA	PP

Model	SHZ-A	SHZ-B
Vacuum port	2	4
Single port flowrate	10L/Min	
Vacuum pressure	0.098Mpa	
Pump Flowrate	60L/Min	
Pump heading	8M	
Pump head material	SS304/PTFE	
Reservoir volume	13L	15L
Motor power	180W	370W
Dimensions	400*280*240mm	
Weight	15Kg	17Kg