

# Laboratory Benchtop Water Circulating Vacuum Pump For Lab Use

Item Number: KV-2



## Introduction

Need a water circulating vacuum pump for your lab or small-scale industry? Our Benchtop Water Circulating Vacuum Pump is perfect for evaporation, distillation, crystallization, and more.

[Learn More](#)

### Detailed material

|         |             |                 |             |           |          |                       |
|---------|-------------|-----------------|-------------|-----------|----------|-----------------------|
| Ejector | check valve | Connecting pipe | three links | pump body | impeller | Air extraction nozzle |
| PP      | Copper+PP   | NR              | PP          | PP        | PA       | PP                    |

|                      |               |       |
|----------------------|---------------|-------|
| Model                | SHZ-A         | SHZ-B |
| Vacuum port          | 2             | 4     |
| Single port flowrate | 10L/Min       |       |
| Vacuum pressure      | 0.098Mpa      |       |
| Pump Flowrate        | 60L/Min       |       |
| Pump heading         | 8M            |       |
| Pump head material   | SS304/PTFE    |       |
| Reservoir volume     | 13L           | 15L   |
| Motor power          | 180W          | 370W  |
| Dimensions           | 400*280*240mm |       |
| Weight               | 15Kg          | 17Kg  |