

Circulating Water Vacuum Pump For Laboratory And Industrial Use

Item Number: KT-VWP



Introduction

Efficient circulating water vacuum pump for labs
- oil-free, corrosion-resistant, quiet operation.
Multiple models available. Get yours now!

[Learn More](#)

Model	SHZ-D(III) Anti-corrosion	SHZ-D(III) PTFE	SHZ-D(III) Stainless Steel	SHZ-D(III) Four Gauges, Four Pumps	SHZ-S(III) Top Water Return	SHZ-95 Vertical Five Pumps
Display Method	Pressure Gauge x2	Pressure Gauge x2	Pressure Gauge x2	Pressure Gauge x4	Pressure Gauge x2	Pressure Gauge x1
Number of Pumping Heads	2	2	2	4	2	5
Max. Vacuum Degree	0.098MPa	0.098MPa	0.098MPa	0.098MPa	0.098MPa	0.098MPa
Single Head Pumping Capacity	10L/min	10L/min	10L/min	10L/min	10L/min	10L/min
Flow Rate	60L/min	60L/min	60L/min	60L/min	90L/min	80L/min
Head	8m	8m	8m	8m	8m	12m
Circulating Water Function	Yes	Yes	Yes	No	Yes	Yes
Pumping Head Material	Copper	PTFE	Copper	Copper	Copper	Copper
Water Tank Volume	15L	15L	15L	15L	15L	35L
Power Supply Voltage	AC 220V 50Hz	AC 220V 50Hz	AC 220V 50Hz	AC 220V 50Hz	AC 220V 50Hz	AC 220V 50Hz
Power	180W	180W	180W	370W	180W	370W
Net Weight	8KG	8KG	9.5KG	13KG	9.5KG	24KG
Product Dimensions	395×295×445	395×295×445	395×295×445	395×295×445	395×295×445	440×340×840