


Custom Ptfе Wafer Holders For Lab And Semiconductor Processing

Item Number: PTFE-38



Introduction

This is a high-purity, custom-machined PTFE (Teflon) holder, expertly designed for the secure handling and processing of delicate substrates like conductive glass, wafers, and optical components.

[Learn More](#)

Holder Size (for Square Sample)	Max Sample Thickness	Recommended Beaker Size
10-15 mm	< 1 mm	250 ml
15-20 mm	< 1 mm	250 ml
21-25 mm	< 1 mm	250 ml
25-30 mm	< 1 mm	250 ml
35-40 mm	< 1 mm	500 ml