

Customizable Pem Electrolysis Cells For Diverse Research Applications

Item Number: KT-PEMT



Introduction

Custom PEM test cell for electrochemical research. Durable, versatile, for fuel cells & CO2 reduction. Fully customizable. Get a quote!

[Learn More](#)

Parameter	Value
Base Material	Metal Plate (Gold-plated or non-gold-plated options available)
Surface Plating	Platinum or non-plated options available
Effective Area	1x1cm, 2x2cm, 5x5cm, 10x10cm (Custom sizes available)
Tubing Diameter	3mm, 6mm (Customizable)
Connector Material	Stainless Steel, PTFE
Custom Processing	Full customization accepted