

# Cylindrical Lab Electric Heating Press Mold

Item Number: PMH



## Introduction

Efficiently prepare samples with Cylindrical Lab Electric Heating Press Mold. Fast heating, high temp & easy operation. Custom sizes available. Perfect for battery, ceramic & biochemical research.

[Learn More](#)

Model	PMH-A	PMH-AG
Heating temp	300°C	500°C
Material	Cr12Mov	
Indenter hardness	HRC60-HRC62	
Sample size	φ10 / φ13 / φ15 / φ20 / φ30 / φ40 mm	
Cavity depth	45 mm	
External dimensions	φ78×138 / 90×138 mm	
Weight(Kg)	2.7	
<b>Other sizes can be customized</b>		