

Cylindrical Press Mold With Scale For Lab

Item Number: PCMC



Introduction

Discover precision with our Cylindrical Press Mold. Ideal for high-pressure applications, it molds various shapes and sizes, ensuring stability and uniformity. Perfect for lab use.

[Learn More](#)

Instrument model	PCMC
Sample shape	
Die heating	Alloy tool steel :Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ10Φ12Φ13Φ15Φ18Φ20 mm (M)
Cavity depth	100mm (N)
Dimensions	Φ53*220mm(L*H)
Weight	4.8Kg
Diagram of hydraulic powder press size	