

Desktop Fast High Pressure Laboratory Autoclave Sterilizer 16L 24L For Lab Use

Item Number: KT-B02



Introduction

The desktop fast steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

Model	KT-T16J	KT-T24J
Sterilization chamber volume	16L $\phi 250 \times 360$ mm	24L $\phi 270 \times 410$ mm
Fixed work stress	0.22Mpa	
Rated working temperature	134°C	
Sterilization temperature selection	105-134°C	
Timing time selection range	0-99min	
Heat uniformity of sterilization room	$\leq \pm 1^\circ\text{C}$	
Power/supply voltage	AC220V.50Hz/2.7KW	
Overall dimensions (mm)	545 \times 328 \times 423	560 \times 350 \times 500
Packing size (mm)	605 \times 390 \times 480	620 \times 410 \times 570
Gross weight/net weight	17/15Kg	20/18Kg