

Desktop Fast Laboratory Autoclave Sterilizer 20L 24L For Lab Use

Item Number: KT-B03



Introduction

The desktop fast steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

Model	KT-XD20D	KT-XD24D
Sterilization chamber volume	20L ϕ 250×420mm	24L ϕ 250×520mm
Fixed work stress	0.22Mpa	
Rated working temperature	134°C	
Sterilization temperature selection	105-134°C	
Timing time selection range	0-99min	
Heat uniformity of sterilization room	≤±1°C	
Power/supply voltage	1.5KW/AC220V.50Hz	
Sterilization tray (mm)	340×200×30(3piece)	400×200×30(3piece)
Overall dimensions (mm)	480×680×455	480×680×455
Packing size (mm)	580×700×500	580×800×560
Gross weight/net weight	43/40Kg	43/40Kg