

Desktop Fast Laboratory Autoclave Sterilizer 35L 50L 90L For Lab Use

Item Number: KT-B04



Introduction

The desktop fast steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items. It efficiently sterilizes surgical instruments, glassware, medicines, and resistant materials, making it suitable for various applications.

[Learn More](#)

| Model | KT-XD35DV | KT-XD50DV | KT-XD90DV |
|---------------------------------------|-------------------|----------------|-----------------|
| Sterilization chamber volume | 35L φ300×500mm | 50L φ340×550mm | 90L φ400×720mm |
| Fixed work stress | 0.22Mpa | | |
| Rated working temperature | 134°C | | |
| Sterilization temperature selection | 105-134°C | | |
| Timing time selection range | 0-99min | | |
| Heat uniformity of sterilization room | ≤±1°C | | |
| Power/supply voltage | AC220V 50Hz/2.7KW | | AC220V 50Hz/4KW |
| Overall dimensions (mm) | 730×550×510 | 770×600×540 | 1000×690×610 |
| Packing size (mm) | 810×680×650 | 880×730×700 | 1070×740×700 |
| Gross weight/net weight | 115/93Kg | 130/106Kg | 200/180Kg |