

# Electric Heating Blast Drying Oven

Item Number: KT-B05



## Introduction

The desktop fast autoclave sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

Model	KT-9040A(S)	KT-9070A(S)	KT-9140A(S)	KT-9240A(S)	KT-9640A(S)
Voltage	AC220V 50Hz				
Temperature range	50°C-300°C				
Temperature fluctuation	±1°C				
Tracking alarm	+5°C				
Power consumption	1000W	1500W	2000W	2500W	4500W
Studio size (mm)	350x350x350	450x350x450	550x450x550	600x500x750	800x800x1000
Dimensions (mm)	500x560x705	650x560x805	750x660x905	800x710x1105	1000x1010x1355