

Laboratory Scientific Electric Heating Blast Drying Oven

Item Number: KT-B05



Introduction

The desktop fast autoclave sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

| | | | | | |
|-------------------------|-------------|-------------|-------------|--------------|----------------|
| Model | KT-9040A(S) | KT-9070A(S) | KT-9140A(S) | KT-9240A(S) | KT-9640A(S) |
| Voltage | AC220V 50Hz | | | | |
| Temperature range | 50°C-300°C | | | | |
| Temperature fluctuation | ±1°C | | | | |
| Tracking alarm | +5°C | | | | |
| Power consumption | 1000W | 1500W | 2000W | 2500W | 4500W |
| Studio size (mm) | 350x350x350 | 450x350x450 | 550x450x550 | 600x500x750 | 800x800x1000 |
| Dimensions (mm) | 500x560x705 | 650x560x805 | 750x660x905 | 800x710x1105 | 1000x1010x1355 |