

Laboratory Diaphragm Vacuum Pump Lab Vacuum Pump

Item Number: LE301



Introduction

Get stable and efficient negative pressure with our Diaphragm Vacuum Pump. Perfect for evaporation, distillation, and more. Low temperature motor, chemical resistant materials, and environmentally friendly. Try it today!

[Learn More](#)

power	parameter
Motor power (W)	180
Power supply	220-240V~, 50Hz
Suction diameter (mm)	∅10×∅6
Maximum vacuum (MPa) /	0.0992/8
Ultimate pressure (mbar)	
Maximum pumping rate (L/min)	25
Dimensions (mm)	310W×225D×168H
Weight (Kg)	10
Operating ambient temperature (°C)	5~35
Ambient relative humidity (%)	≤ 70
Enclosure protection class	IP20
Anti-pollution	level 2

Inlet and outlet nozzle connection hose specifications: (outer diameter × inner diameter) (mm) Rubber hose: ∅12×∅6