

Double Plate Heating Press Mold For Lab

Item Number: PMD



Introduction

Discover precision in heating with our Double Plate Heating Mold, featuring high-quality steel and uniform temperature control for efficient lab processes. Ideal for various thermal applications.

[Learn More](#)

Instrument model	PMD
Press the shape of the sample	
Heating temperature	Room temperature-300C
Mold material	Alloy tool steel:Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	
Sample thickness	0.02-0.1mm(N)
External dimension 1	180*180*130mm(L*W*H)
External dimension 2	200*200*130mm(L*XW*H)
External dimension 3	300*300*130mm(L*W*H)
Mold weight	32Kg38Kg88Kg
Dimensional diagram	