

Lab Scale Rotary Single Punch Tablet Press Machine Tdp Tablet Punching Machine

Item Number: TPM-01



Introduction

This machine is a single-pressure automatic rotating, continuous tableting machine that compresses granular raw materials into various tablets. It is mainly used for tablet production in the pharmaceutical industry, and is also suitable for chemical, food, electronics and other industrial sectors.

[Learn More](#)

| | |
|----------------------|----------------|
| Die quantity | 5 pcs |
| Pressure | 60KN |
| Max. Tablet diameter | 18mm |
| Max. Tablet height | 6mm |
| Max. Filling height | 15mm |
| Disc rotate speed | 30rpm |
| Capacity | 9000pcs/h |
| Motor power | 2.2Kw |
| Dimension | 960*660*1230mm |
| Weight | 320kg |

1. Excess powders can be recycled by a vacuum cleaner(Optional)
2. Tablet die can customized according to specific requirement