

# Lab Scale Rotary Tablet Punching Machine

Item Number: TPM-01



## Introduction

This machine is a single-pressure automatic rotating, continuous tableting machine that compresses granular raw materials into various tablets. It is mainly used for tablet production in the pharmaceutical industry, and is also suitable for chemical, food, electronics and other industrial sectors.

[Learn More](#)

Die quantity	5 pcs
Pressure	60KN
Max. Tablet diameter	18mm
Max. Tablet height	6mm
Max. Filling height	15mm
Disc rotate speed	30rpm
Capacity	9000pcs/h
Motor power	2.2Kw
Dimension	960*660*1230mm
Weight	320kg

1. Excess powders can be recycled by a vacuum cleaner(Optional)
2. Tablet die can customized according to specific requirement