

# Electrolytic Electrochemical Cell With Five-Port

Item Number: ELC



## Introduction

Streamline your laboratory consumables with Kintek's Electrolytic Cell with five-port design. Choose from sealed and non-sealed options with customizable electrodes. Order now.

[Learn More](#)

Specification	10ml~1000ml
Applicable temperature range	0~60°C
Sealed form	Thread / Apron
Material	Boron glass, PTFE
Openings in the electrolytic cell	Three electrode holes (6mm) and Two air holes (3mm), can be customized
Specification	10ml~1000ml
Applicable temperature range	0~60°C
Material	Boron glass, PTFE
Openings in the electrolytic cell	Three electrode holes (6mm), can be customized