

Flat Corrosion Electrolytic Electrochemical Cell

Item Number: ELEFC



Introduction

Discover our flat corrosion electrolytic cell for electrochemical experiments. With exceptional corrosion resistance and complete specifications, our cell guarantees optimal performance. Our high-quality materials and good sealing ensure a safe and durable product, and customization options are available.

[Learn More](#)

Specifications	350ml, can be customized
Applicable temperature range	0 ~ 70°C
Sealing form	TSilicone rubber gasket
Material	boron glass + PTFE
Hole	three grinding mouths + two inner circulation pagoda mouths
Specifications	350ml, can be customized
Applicable temperature range	0 ~ 100°C
Material	boron glass + PTFE
Hole	Three grinding mouth two circulation + water bath