

# Electrolytic Electrochemical Cell Gas Diffusion Liquid Flow Reaction Cell

Item Number: ELCG



## Introduction

Looking for a high-quality gas diffusion electrolysis cell? Our liquid flow reaction cell boasts exceptional corrosion resistance and complete specifications, with customizable options available to suit your needs. Contact us today!

[Learn More](#)

Airway type	snake-shaped airway / back-shaped airway / special-shaped custom
Cell material	optional PTFE / PEEK / PP / plexiglass / nylon