

Glass Oscillating Beads Laboratory Zeolite Transparent Glass Balls

Item Number: KTG-4



Introduction

Glass vibration beads, commonly used in laboratory settings, are transparent glass balls designed to prevent zeolite formation.

[Learn More](#)

Diameter	Whole package weight	Packaging
2.5mm	18.68g	1000 pcs/pack
3mm	39.92g	1000 pcs/pack
4mm	81.58g	1000 pcs/pack
4.5mm	123.57g	1000 pcs/pack
5mm	159.39g	1000 pcs/pack
6mm	281.06g	1000 pcs/pack
7mm	451.38g	1000 pcs/pack